

MVE CryoCarts are designed for the loading of biological samples into canes, boxes, racks or frames. When used as a portable workbench, the unit will provide a safe and controlled environment for your samples for up to eight hours with the lid open. MVE CryoCart is ideal for transporting large quantities of samples from one tank to another within the same facility.

## Features Include:

- Approximate hold time of 18 hours with lid on
- Omega temperature monitor with type T Thermocouple
- Traceable calibration certificate with RS232 PC interface
- LN2 transfer hose connection with ball valve and relief valve
- Available without plumbing and temperature recorder.

## Product Specifications

### UNIT CAPACITIES

LN2 Capacity To Top of Platform *L* 53

### UNIT DIMENSIONS

<b>Overall Length in (mm)</b>	51.7 (1313)
<b>Overall Width in (mm)</b>	17.6 (447)
<b>Overall Height in (mm)</b>	38.9 (988)
<b>Lift Overall Height in (mm)</b>	36.3 (922)
<b>Inside Length in (mm)</b>	41.0 (1041)
<b>Inside Width in (mm)</b>	13.2 (335)
<b>Inside Depth in (mm)</b>	15.4 (391)
<b>Platform Height in (mm)</b>	6.0 (152)
<b>Footprint in (mm)</b>	17.6x51.7 (447x1313)
<b>Weight Empty lb. (kg)</b>	241 (109)
<b>Weight w/LN2 Filled To</b>	
<b>Top Of Platform lb. (kg)</b>	356 (161)
<b>Approximate Hold Time With Lid On hrs</b>	18
<b>Number of Racks (15-2)</b>	2



## ONE (1) Year Standard Warranty

**Chart BioMedical**  
 2200 Airport Industrial Dr., Ste. 500  
 Ball Ground, GA 30107  
 Ph 770-721-7759 • Toll Free 1-800-482-2473  
 Fax 770-721-7758  
 customerservice.usa@chartindustries.com



www.chartbiomed.com

**Chart BioMedical Ltd.**  
 Ph +44(0) 1344 403100 • Fax +44(0) 1344 429224  
 customerservice.europe@chartindustries.com  
**Chart Japan Co., Ltd.**  
 Ph + (03) 5776 2670 • Fax + (03) 5776 2676  
 customerservice.japan@chartindustries.com  
**Chart Australia Pty Ltd.**  
 Ph + 61 2 974 94 333 • Fax +61 2 974 94 666  
 customerservice.australia@chartindustries.com